

TRISTAN TECHNOLOGIES, INC

---

Magnetic and Cryogenic Device solutions

Bio Systems  
Information Brochure

TRISTAN TECHNOLOGIES, INC.

# ***Bio Systems Information Brochure***

---

© Tristan Technologies, Inc  
6191 Cornerstone Court East • Suite 107  
San Diego, CA 92121  
Phone 858.550.2700 • Fax 858.550.2799

# ***Table of Contents***

SYSTEM INFORMATION AND ADVANTAGES.....	2
BabySQUID NEONATAL BIOMAGNETOMETER .....	5
ARTEMIS123®.....	7
A) EchofMCG™ SQUID SYSTEM FOR FETAL CARDIAC MEASUREMENTS.....	10
B) TruckSQUID™ SYSTEM FOR MOBILE FETAL MAGNETOCARDIOGRAPHY (fMCG).....	11
BIOMAGNETIC LIVER SUSCEPTOMETER.....	12
INTESTINAL ISCHEMIA SYSTEM.....	14
SpineSQUID™ SPINAL CORD MEASUREMENT SYSTEM...	17
MODEL 601 SINGLE CHANNEL GRADIOMETER SYSTEM..	18
INVERTED SQUID MICROSCOPE.....	19
INFANT MEG SYSTEM- MagView™ .....	20

# ***System information and advantages***

## ***Biomagnetic measurements provide a number of advantages compared to electrical measurements:***

- Biomagnetism is non-invasive. The detection system does not contact the subject. The non-invasive nature of biomagnetism makes it an inherently safe procedure and minimizes subject preparation time.
- Insulating barriers such as the skull, varying layers of tissue, anatomical open spaces, do not attenuate or distort magnetic fields. Electrical signals are distorted by the varying resistive layers between the signal source and the surface skin.
- SQUID magnetometers will measure the vector component(s) of the magnetic field. Thus localization is much easier than with electrical measurements, which only measure scalar voltages.
- Magnetic measurements can be made for which there are no electrical analogs. These include measurements of static magnetic fields, measurements of the magnetic susceptibility and measurements where an invasive procedure is not possible (e.g., fetal cardiography).
- Because of the superconducting nature of SQUID measurements, true dc response and flat phase response are available.

## ***Instrumentation***

The strength of biomagnetic signals is many orders of magnitude smaller than even the earth's magnetic field, which is 1/2 Gauss or 50 microtesla. The signal strengths associated with biomagnetism (Fig. 1) require the use of extremely sensitive detection systems. The units in this figure are femtotesla, 1 fT = 10<sup>-15</sup> tesla. The only instrument with the required sensitivity and bandwidth is the SQUID magnetometer.

The components of a SQUID magnetometer (Fig. 2) typically consist of the following: a detection coil, which senses changes in the external magnetic field and transforms them into an electrical current; an input coil which transforms the resulting current into a magnetic flux in the SQUID sensor; electronics which transform the applied flux into a room temperature voltage output; and acquisition hardware and software for acquiring, storing and analyzing data. Both the SQUID amplifier and the detection coils are superconducting devices. Thus some type of refrigerant (liquid helium or liquid nitrogen) or refrigeration device (cryocooler) is needed to maintain the SQUID and detection coil in the superconducting state. Additional signal conditioning electronics may be needed to improve signal-to-noise.




Figure 1: Typical amplitudes and frequencies of various biomagnetic signals. shaded regions indicate the limits of currently available SQUID sensors.




Figure 2: Block diagram of a SQUID magnetometer

The SQUID sensor and electronics package can be considered as a black box that acts like a magnetic field-to-voltage converter and amplifier with extremely high gain. In addition, it offers extremely low noise, high dynamic range, excellent linearity, flat phase response and a bandwidth that can extend from dc to beyond 100 kHz, capabilities that no other single sensor offers.

The type of SQUID sensor and detection coil configuration is dependent on what is to be measured. Figure 1 also shows the capability of both low temperature (requiring liquid helium temperatures, and referred to as LTS) and high temperature (requiring liquid nitrogen temperatures, and referred to as HTS) SQUID magnetometers. Tristan biomagnetic measurement systems make use of either Tristan's LSQ/20 LTS dc SQUID sensor or the HTM-8 HTS dc SQUID sensor. The input coil for an LTS SQUID is normally fabricated from flexible superconducting NbTi wire. The inherent anisotropic nature of HTS SQUIDs requires that the input coils be planar. Typically HTS magnetometers are available only as pure magnetometers.

Another factor to be considered is the detection coil configuration. Conceptually, the easiest input circuit to consider for detecting changes in magnetic fields is a pure magnetometer (Fig. 2). However, magnetometers are extremely sensitive to all magnetic signals in the environment. This may be acceptable if one is measuring ambient fields. However, if the magnetic signal of interest is weak, then environmental magnetic interference may prevent measurements. If the signal source is close to the detection coil, then a gradiometer coil may allow a weak signal to be measured. Figure 3 shows the relative noise rejection for 1st and 2nd derivative gradiometers. The figure insert shows a first order gradiometer, consisting of two coils connected in series but wound in opposite senses, and separated by a distance "b", called the gradiometer baseline. A uniform magnetic field (e.g., from a distant environmental source) would couple equal but opposite quantities of flux into the two coils, resulting in zero net flux in the gradiometer, or zero signal. However, signal sources that are close to the lower coil (relative to the baseline, or separation

between coils) would couple significantly more flux into the lower coil than into the upper coil; this would result in a net flux in the gradiometer and hence the signal would be detected.





Figure 3: Response of gradient coils relative to magnetometer response (1/r<sup>3</sup> suppressed)

For objects that are close (relative to the gradiometer baseline), the gradiometer acts as a pure magnetometer, while rejecting more than 99% of the magnetic signals coming from distant objects. In essence, the gradiometer acts as a “compensated” magnetometer.

Normally, SQUID magnetometers (and gradiometers) map the axial (B<sub>Z</sub>) component of the magnetic field. Obviously, using three sensors, it is possible to monitor all three vector components of the magnetic field. Additional channels of SQUID sensors can be used to provide reference channels for electronic balancing. Portions of the reference magnetometer responses are summed electronically with the detection coil(s) output to reject common mode signals from distant noise sources. Electronic balancing can be used to create an HTS axial gradiometer from two HTS magnetometers.

# 1. **BabySQUID NEONATAL BIOMAGNETOMETER**

The BabySQUID Neonatal Biomagnetometer is a new investigational tool for neurological impairments of pre-term and full-term infants. babySQUID® (pdf) measures and maps brain activity non-invasively at the



bedside.

## **CLINICAL ADVANTAGES of babySQUID®**

- Map the sites and dynamics of sensory functions
- Map seizure and inter-ictal activity for epilepsy monitoring
- Assay stages of nervous system development
- Monitor recovery from trauma
- Detect effects of earlier hypoxic and intracranial injury

**More newborns survive even with neurological disabilities.**

<b>CONDITION</b>	<b>INCIDENCE</b>
perinatal asphyxia	between 2 – 47 per1000
hypoxemic-ischemic encephalopathy	between 3 – 8 per 1000
moderate-to-severe cerebral palsy (post neonatal)	between 1 – 3 per 1000
periventricular white matter injury	240 per 1000 (for gestational age < 38 weeks)

## **TECHNICAL ADVANTAGES OF BabySQUID®**

- Unprecedented spatial resolution and sensitivity.
- A dense array of closely spaced sensors is located just below the outer surface of a headrest.
- The sensor noise is  $< 10 \text{ fT}/\sqrt{\text{Hz}}$  for the detection coils.
- BabySQUID® has an order of magnitude better sensitivity to neuronal sources than conventional whole-head MEG systems
  - o Sensitive enough to measure spontaneous neuronal activity and evoked activity of the cortex of the newborns in real time without signal averaging
  - o Spatial resolution four times greater than existing whole-head MEG sensors
- In comparison, EEG signals are significantly distorted by skull defects (fontanelles and sutures) unique to the human neonates. These skull defects can obscure the asymmetry of the signals, especially when the generator is deep, making it difficult to determine the location of the epileptiform tissue when it cannot be easily visualized by CT or MRI.

## **PRINCIPLES OF PRODUCTION**

Superconducting amplifiers (SQUIDs) are used to amplify magnetic signals detected by a large array of small detection coils. The detection apparatus is kept at cryogenic temperatures by a vacuum insulated vessel (dewar). Tristram developed fabrication methods allow the detection coils to be placed extremely close to the patient without loss of sensitivity or risk to the patient. Additional information on neuromagnetic instrumentation can be found at [Chapt2.pdf](#). The babySQUID® takes advantage of the fact that the infant's scalp and skull are thin. This makes it possible to measure MEG signals at a distance of only about 5-6 mm from the brain surface. This shorter distance results in a significant increase in amplitude of MEG signals from the newborns, since the magnetic field is inversely proportional to the square of the distance. The shorter distance and the high density of detectors also results in higher spatial resolution.




## **2 ARTEMIS123®**



### ***SYSTEM DESCRIPTION, PRINCIPLES OF OPERATION***


Like adult Magnetoencephalography (MEG) systems, Artemis123® uses superconducting sensors to non-invasively detect and map magnetic fields generated by cortical neural activity. However, Artemis123® takes advantage of the fact that the infant's scalp and skull are very thin. Tristan's fabrication methods put the sensing coils very close to the infant brain's sources of activity, even though SQUIDs must operate in an ultra-cold liquid helium environment. The net result is a significant increase in amplitude of neonate MEG signals. Also, the high density of detectors results in higher spatial resolution compared to adult whole-head MEG.



Mapping of sites and dynamics of sensory functions – auditory, somatosensory, and visual modalities.

## **SYSTEM COMPONENTS**

- Sensor/Cradle/Bed on mobile cart – easily accessed height
- Power supplies and computer on companion mobile cart to minimize noise
- Subject Tracking – optical tracking system updates movement at 30 Hz with ½ mm accuracy
- Part-wise mapping or optional optical one-click 3D imaging system
- Assay stages of nervous system development



Somatic evoked magnetic field (SEF) obtained from a 7-month old as a function of number of averages from N=4 to 173 epochs. The waveforms are the differences of the SEF at two field extrema. This shows that a small number of averages are needed to acquire SEF data. (data acquired using a Tristan babySQUID® system).

## **UNIQUE FEATURES OF ARTEMIS123®**

- Superior spatial resolution and sensitivity
- Significantly more sensitive to neuronal sources than conventional whole-head MEG systems
- Similar or better spatial resolution compared to existing whole-head MEG sensors
- Better spatial resolution than EEG (EEG signals are distorted by skull defects (fontanels and sutures), making it difficult to localize epileptiform tissue)

- Rapid scanning: a typical clinical scan can be completed within thirty minutes
- Anti-vibration construction; infant motion will not cause vibrational artifacts
- Sensor noise level  $< 10 \text{ fT/ Hz}$
- A dense array of closely-spaced sensors located just below the outer surface of a headrest.
- Allows simultaneous measurement of the occipi-tal area and parietal and temporal areas
- Includes position tracking device and software, permitting measurements during sleep or relatively quiescent wakefulness
- The measurement cradle and companion electronics cart are portable and can be wheeled in and out of elevators, obstetric suites and neo-nate ICUs

### 3. A) **EchofMCG™ SQUID SYSTEM FOR FETAL CARDIAC MEASUREMENTS**



The **EchofMCG™** is a unique multi-channel vector fMCG system integrated with an echocardiography system, capable of simultaneously performing fetal magnetocardiography (fMCG) and echocardiography (echo/Doppler). Such linking of echo/Doppler and fMCG would allow the clinician to analyze the fetus rapidly for both hemodynamic as well as electrophysiologic abnormalities, such as fetal arrhythmias. This will be the world first clinical modality to provide full characterization of the intrauterine condition of fetuses with life threatening heart conditions.

**SQUID = Superconducting Quantum Interference Device.**

- Simultaneously detection of both ultrasound and magnetic fields associated with fetal cardiac electrical activity
- Vector gradiometer design to maximize captured fMCG information
- Small probe profile for easy positioning in close proximity to fetus
- Gantry movement offers four degrees of freedom for patient accommodation

#### 21-Channel Signal Array

- 7 vector gradiometer elements ( $dB_x/dz$ ,  $dB_y/dz$ ,  $dB_z/dz$ ) arranged in a hexagonal array,
- with ~8 cm baselines to allow for variability in the source depth
- Coil-coil spacing of ~4 cm and an overall tail diameter of ~10 cm



#### The EchofMCG™ System Shown With

- 7 Ultrasound probe mounted on the fMCG dewar probe

#### **EchofMCG™ SYSTEM ADVANTAGES**

- Vector field mapping capability.
- Ability for deep source detection
- Synchronized fMCG and fetal ultrasound measurement.
- Advanced data processing based on spatial filtering and ICA.

## ***B) TruckSQUID™ SYSTEM FOR MOBILE FETAL MAGNETOCARDIOGRAPHY (fMCG)***


The TruckSQUID™ is a unique system for fetal magnetocardiography (fMCG) measurements. It allows the clinician to analyze a fetus rapidly for electrophysiologic abnormalities such as fetal arrhythmias. This system is the first clinical mobile system that provides full intrauterine characterization of a fetus with life-threatening heart conditions.



### ***TruckSQUID™ SYSTEM CHARACTERISTICS:***

- Vector field mapping capability
- Deep source detection capability Liquid Helium
- Dewar hold time 5-7 days Windows-based acquisition and display software
- Advanced data processing based on spatial filtering
- ICA optional Ultrasound probe

## 4. BIOMAGNETIC LIVER SUSCEPTOMETER



The BLS Liver Iron Stores Measurement System is designed for measuring fields from paramagnetic materials in the body, such as hepatic iron stores in the liver. Measurements are made by determining the change in magnetic field at the detector as the subject is moved into and away from the sensitive region of the detector. A small magnetic field is applied during these measurements by a self-contained superconducting magnet. To simulate the presence of the body during the measurements, water approximating the natural diamagnetism of the body is located between the sensor and the body.

The system includes dual channel axial gradiometers (3rd channel optional), superconducting magnet and power supply, a liquid helium dewar and gantry, water bag and reservoir, movable bed, a data acquisition and analysis system, and all necessary accessories. As with all Tristan systems, an on-site training course in the proper use of the system is available.

### APPLICATIONS

The most relevant applications of Biomagnetic Liver Susceptometry (BLS) are related to iron overload diseases such as hereditary hemochromatosis and siderosis caused by blood transfusions. To date, the following applications have been demonstrated:

- Monitoring iron overload in patients with transfusional siderosis (genetic  $\beta$ -thalassemia major and sickle cell disease, or other transfusion dependent anemias) for the onset or intensification of chelation therapy and during this therapy.
- Assessment of iron overload in patients scheduled for Interferon alfa therapy in viral liver infections such as Hepatitis B or C.
- Assessing iron overload in patients with  $\beta$ -thalassemia scheduled for bone marrow transplantation (BMT) or monitoring iron overload after BMT during iron depletion therapy.
- Assessment of the long-term efficacy of different iron chelators under study.
- Diagnosis of hereditary hemochromatosis and assessment of the degree of iron overload in known hereditary hemochromatosis.
- Monitoring liver iron concentration in the initial assessment and long term phlebotomy therapy of hereditary hemochromatosis.

## **METHODOLOGY**

Non-invasive Biomagnetic Liver Susceptometry (BLS) exploits the effects of magnetism and superconductivity. Biological materials such as ferritin and hemosiderin are weakly attracted to an applied magnetic field (paramagnetic behavior) while water and body tissue are very weakly repelled (diamagnetic).

Ferromagnetic materials e.g., nickel and steel, are strongly attracted to applied fields. No naturally occurring human tissue is ferromagnetic. In the BLS method, a weak magnetic field of 0 - 20 millitesla is generated within the body tissue by an external superconducting field magnet, similar to that used in a MRI scanner, but a hundred times weaker. The applied fields are measured by a superconducting magnetometer known as a SQUID (Superconducting Quantum Interference Device). The SQUID sensing system has the ability to measure distortions in the magnetic field at the part per billion levels.

When an organ, such as the liver, is placed in a magnetic field, it will slightly distort the applied field. If the liver is normal or anemic, the local field will be reduced slightly. If the liver is iron overloaded, the local field will be enhanced. Hence the change in the detected magnetic field is directly related to the iron concentration in the liver. To minimize the body's contribution to the distortion in magnetic field, a small bag of water is placed between the detector and skin surface. Since the susceptibility of body tissue is close to that of water, the resultant measurement is essentially that of a magnetized liver (or spleen) moving in a magnetic field within a uniform (diamagnetic) environment; the only change seen by the detection coils is due to the liver (or spleen) itself.


For higher accuracy, our software removes the actual contribution of overlying tissues (skin, bone, muscle, fat, etc.). This gives the iron concentration of the liver (or spleen) alone, allowing accurate measurements for obese patients and normal patients with atypical liver/spleen depths.

To date, the BLS method has been applied to organs such as livers and enlarged spleens (> 300 ml) with a total error of  $[Fe] = 0.05 - 0.4 \text{ mg/g tissue (wet weight)}$ . Repeatability (serial measurements over three weeks) on single subjects of better than 95 % has been demonstrated.



## 5. **INTESTINAL ISCHEMIA SYSTEM**

Tristan Technologies fabricates a high sensitivity, multi-channel SQUID magnetometer system for measuring electromagnetic activity in the human intestine. Presently, intestinal ischemia is difficult to diagnose, and is usually fatal. SQUID sensors can detect the magnetic fields produced by the BER (basic electrical rhythm) of the human intestine. The frequency of the BER signals changes under ischemia — the frequency of BER intestinal signals are ~10 cpm (cycles per minute).



Magnetic measurements provide improved signal-to-noise over the currently more typical cutaneous electrode measurements of electric potential. In contrast to the measurements of voltages on the skin surface, magnetic signals are not attenuated or redirected by the multiple layers of varying electrical resistivity tissues separating the intestine from the skin surface. With multi-channel magnetic measurements, vector projection analysis techniques allow focusing on the signals of interest, distinguishing them from the many other biomagnetic and environmental signals present. Other less serious intestinal disorders, such as Crohn's disease, ulcerative colitis, and irritable bowel, are also difficult to diagnose; their diagnoses may be improved with this system.

- Non-Invasive — no contact between instrument and abdominal wall.
- Magnetic measurements superior to electric
- Signals not attenuated or redirected by the multiple layers of tissue separating skin from intestine
- Improved signal-to-noise.
- Detect signal changes before pathological damage.
- Useful information in short time periods — extensive patient preparation or analysis not required.







Figure 1. Schematic drawing: Model 637 Intestinal Ischemia System.

## ***ELEMENTS IN THE MODEL 637 INTESTINAL ISCHEMIA SYSTEM***

- 29 magnetic field sensing channels, < 20 mm from sensor surface, distributed over
- Large (296 cm<sup>2</sup>) area of coverage
- Or, intermediate (82 cm<sup>2</sup>) area of coverage (set at Tristan facility)
- 8 magnetic sensing channels, in a tensor array, monitoring environmental magnetic noise.



This system is a multi-channel system with 29 detection coils (19 axial coils, 10 vector coils) designed to measure the Basic Electric Rhythm (BER) associated with intestinal activity. The specific application is detection of Mesenteric Ischemia, a life threatening condition with no conventional reliable method of diagnosis. Pre-clinical trials in partnership with Vanderbilt University are underway. The system features coil-in-vacuum construction and the unique ability to vary the position of the detection coils. This allows the researcher to adjust the spatial frequency measurement capability of the system. Like the model 619, the 'model 637' includes an 8-element tensor noise reduction scheme. Designed to operate in a clinical setting, the model 637 operates in an unshielded environment. Tristan's experience with coil-in-vacuum design is critical for sensors that are both portable and adaptable to measurement at varying orientations.

## 6. **SpineSQUID™ SPINAL CORD MEASUREMENT SYSTEM**

Tristan has built and delivered a fully integrated 63 channel magnetic source imaging system for non-invasive measurements of spinal cord activity and source localization. The system is adaptable for humans or animals. Because spinal signals are action potentials, the system is designed to acquire data in excess of 100,000 samples per second on each of its 80 channels (including reference channels), more than an order of magnitude faster than conventional MEG devices. For this project, Tristan devised a novel high speed data acquisition and monitoring system capable of acquiring and storing more than 10 minutes of continuous spine data, and simultaneously retrieving and reviewing a data set collected previously.



### **ADVANTAGES**

- Tailor-made complex shape Dewars Sophisticated cryogenic dewar construction
- Noise reduction software and hardware
- Orthogonally oriented vector ( $B_x$ ,  $B_y$ ,  $B_z$ ) detection coils Asymmetric gradiometers for improved sensitivity Multi-axis dewar gantries
- Magnetically quiet motorized patient beds Vector and tensor reference arrays for noise cancellation Software for data analysis in LabView™ including FFT, digital filtering and real-time review
- Software compatibility with standard source analysis packages (BESA, EMSE) and with MATLAB.
- Integrated hardware and software for positioning and tracking the subject: including 3D optical positioning camera system.

## **7. MODEL 601 SINGLE CHANNEL GRADIOMETER SYSTEM**

The 601 is a single channel LTS (liquid helium) SQUID gradiometer system. Its components consist of a Cryogenic Probe with liquid helium level sensor, a 1st order axial (dBz/dz) detection coil, iMAG® LTS SQUID and electronics (1 channel) and a Model BMD-6 Liquid Helium Dewar. With a 1 cm detection coil, sensitivities approaching 10 fT/ÖHz are possible. The BMD-6 dewar allows the detection coils to be placed within 10 mm of room temperature. System components:

- 1st order axial detection coil, nominal 1 cm diameter, 2% balance
- Cryogenic Probe with liquid helium level sensor
- Model LSQ/20 LTS dc SQUID
- Model BMD-6 Liquid Helium Dewar
- Model iMC-303 Cryogenic Control Unit
- Model iFL-301-L Flux-Locked Loop
- Model CC-6 six meter fiber-optic composite connecting cable
- Manual and accessory pack

### **SAMPLE CONFIGURATION**

- 2mm diameter 1<sup>st</sup> order gradiometers with 1 cm baseline.
- 2mm coil to coil separation.
- 2mm offset from room temperature outer dewar surface.



## 8. ***INVERTED SQUID MICROSCOPE***



Tristan Technologies has developed a prototype Inverted SQUID (Superconducting Quantum Interference Device) Microscope for neuroscience research. The target signal levels are much weaker (100-500 fT, fT = 10-15 Tesla) than signals in the area of non-destructive evaluation (> 1 pT, pT= 10-12 Tesla) where SQUID microscopes have been used previously. The term “inverted” is adopted because the microscope is similar to an inverted optical microscope except the objective lens is replaced by an array of superconducting miniature magnetic field sensing coils.

The microscope is useful for other applications that include measurements of:

- (1) Electrical currents from single neurons and glial cells in culture
- (2) Efficiency of bonding of antigens and magnetically tagged antibodies (immunoassay)
- (3) Movements and conformational changes of a small number of magnetically tagged molecules in a cell for studying signaling pathways.

The inverted SQUID microscope is useful in both academic setting and industry for understanding the electrophysiology of small cells that are difficult to study with electrodes, for drug discovery and for studying second-messenger systems.

## 9. **INFANT MEG SYSTEM- MagView™**



The Tristan MAGView™ Biomagnetometer features whole head coverage for a helmet designed to fit a 50 cm circumference head. It is used to non-invasively measure weak magnetic fields produced by electrical activity from the brain of infants and children. The system consists of the following principal components: the sensor, a mobile patient bed, an electronics cart containing SQUID electronics, an external electronics rack for power supplies and data acquisition hardware, and a computer. The patient bed, sensor, and SQUID electronics rack are designed to fit inside a magnetically shielded room (MSR).

### **ADVANTAGES**

- Superior spatial resolution and sensitivity
- Significantly more sensitive to neuronal sources than conventional whole-head MEG systems
- Similar or better spatial resolution compared to existing whole-head MEG sensors
- Better spatial resolution than EEG. EEG signals are distorted by skull defects (fontanels and sutures), making it difficult to localize epileptiform tissue
- Rapid scanning: a typical clinical scan can be completed within thirty minutes
- Anti-vibration construction; infant motion will not cause vibrational artifacts
- Sensor noise level  $< 10 \text{ fT}/\sqrt{\text{Hz}}$
- A dense array of closely-spaced sensors located just below the outer surface of a helmet.
- Allows simultaneous measurement of the occipital area, parietal areas, and temporal areas
- Includes position tracking device and software, permitting measurements during sleep or relatively quiescent wakefulness




## **SYSTEM DESCRIPTION**

Like adult Magnetoencephalography (MEG) systems, MAGView™ uses superconducting sensors to non-invasively detect and map magnetic fields generated by cortical neural activity. However, MAGView™ takes advantage of the fact that the infant's scalp and skull are very thin. Tristan's fabrication methods put the sensing coils very close to the infant brain's sources of activity, even though SQUIDs must operate in an ultra-cold liquid helium environment. The net result is a significant increase in amplitude of neonate MEG signals. Also, the high density of detectors results in higher spatial resolution compared to adult whole-head MEG.

The MAGView™ signal detector channels are specified to have a noise level and sensitivity to magnetic fields of at least  $10 \text{ fT}/\sqrt{\text{Hz}}$  or better on average. Ambient magnetic fields in a typical hospital environment are generally much greater than this sensitivity, and in many cases, the system will be operated within a magnetically shielded room to enable measurements with the full sensitivity capability.

## **SQUID SENSOR ARRAY**



ARRAY ARRANGEMENT WITHIN THE DEWAR

- 200 to 400 channel sensors within the helmet
- Magnetometer detectors
- Reference channels for ambient noise reduction
- Coil-in-vacuum configuration for superconducting coil array and SQUIDs
- Coil-to-surface gap  $\sim 6 \text{ mm}$
- Average system white noise  $< 10 \text{ fT}/\sqrt{\text{Hz}}$  in magnetically quiet environment
- Helmet designed for whole head coverage, with 50 cm circumference
- Helmet positioned at a height between 30-36" from MSR floor
- Subjects measured in a supine position

# Non-Bio Systems Information Brochure



TRISTAN TECHNOLOGIES, INC.

# ***Non-Bio Systems Information Brochure***

---

© Tristan Technologies, Inc  
6191 Cornerstone Court East • Suite 107  
San Diego, CA 92121  
Phone 858.550.2700 • Fax 858.550.2799

# ***Table of Contents***

SYSTEM INFORMATION AND ADVANTAGES.....	2
LABORATORY APPLICATIONS.....	5
GEOPHYSICAL APPLICATIONS.....	5
NON-DESTRUCTIVE TEST & EVALUATION.....	5
MEDICAL APPLICATIONS.....	6

# ***System information and advantages***

## ***Biomagnetic measurements provide a number of advantages compared to electrical measurements:***

- Biomagnetism is non-invasive. The detection system does not contact the subject. The non-invasive nature of biomagnetism makes it an inherently safe procedure and minimizes subject preparation time.
- Insulating barriers such as the skull, varying layers of tissue, anatomical open spaces, do not attenuate or distort magnetic fields. Electrical signals are distorted by the varying resistive layers between the signal source and the surface skin.
- SQUID magnetometers will measure the vector component(s) of the magnetic field. Thus localization is much easier than with electrical measurements, which only measure scalar voltages.
- Magnetic measurements can be made for which there are no electrical analogs. These include measurements of static magnetic fields, measurements of the magnetic susceptibility and measurements where an invasive procedure is not possible (e.g., fetal cardiography).
- Because of the superconducting nature of SQUID measurements, true dc response and flat phase response are available.

## ***Instrumentation***

The strength of biomagnetic signals is many orders of magnitude smaller than even the earth's magnetic field, which is 1/2 Gauss or 50 microtesla. The signal strengths associated with biomagnetism (Fig. 1) require the use of extremely sensitive detection systems. The units in this figure are femtotesla, 1 fT = 10<sup>-15</sup> tesla. The only instrument with the required sensitivity and bandwidth is the SQUID magnetometer.

The components of a SQUID magnetometer (Fig. 2) typically consist of the following: a detection coil, which senses changes in the external magnetic field and transforms them into an electrical current; an input coil which transforms the resulting current into a magnetic flux in the SQUID sensor; electronics which transform the applied flux into a room temperature voltage output; and acquisition hardware and software for acquiring, storing and analyzing data. Both the SQUID amplifier and the detection coils are superconducting devices. Thus some type of refrigerant (liquid helium or liquid nitrogen) or refrigeration device (cryocooler) is needed to maintain the SQUID and detection coil in the superconducting state. Additional signal conditioning electronics may be needed to improve signal-to-noise.




Figure 1: Typical amplitudes and frequencies of various biomagnetic signals. shaded regions indicate the limits of currently available SQUID sensors.




Figure 2: Block diagram of a SQUID magnetometer

The SQUID sensor and electronics package can be considered as a black box that acts like a magnetic field-to-voltage converter and amplifier with extremely high gain. In addition, it offers extremely low noise, high dynamic range, excellent linearity, flat phase response and a bandwidth that can extend from dc to beyond 100 kHz, capabilities that no other single sensor offers.

The type of SQUID sensor and detection coil configuration is dependent on what is to be measured. Figure 1 also shows the capability of both low temperature (requiring liquid helium temperatures, and referred to as LTS) and high temperature (requiring liquid nitrogen temperatures, and referred to as HTS) SQUID magnetometers. Tristan biomagnetic measurement systems make use of either Tristan's LSQ/20 LTS dc SQUID sensor or the HTM-8 HTS dc SQUID sensor. The input coil for an LTS SQUID is normally fabricated from flexible superconducting NbTi wire. The inherent anisotropic nature of HTS SQUIDS requires that the input coils be planar. Typically HTS magnetometers are available only as pure magnetometers.

Another factor to be considered is the detection coil configuration. Conceptually, the easiest input circuit to consider for detecting changes in magnetic fields is a pure magnetometer (Fig. 2). However, magnetometers are extremely sensitive to all magnetic signals in the environment. This may be acceptable if one is measuring ambient fields. However, if the magnetic signal of interest is weak, then environmental magnetic interference may prevent measurements. If the signal source is close to the detection coil, then a gradiometer coil may allow a weak signal to be measured. Figure 3 shows the relative noise rejection for 1st and 2nd derivative gradiometers. The figure insert shows a first order gradiometer, consisting of two coils connected in series but wound in opposite senses, and separated by a distance "b", called the gradiometer baseline. A uniform magnetic field (e.g., from a distant environmental source) would couple equal but opposite quantities of flux into the two coils, resulting in zero net flux in the gradiometer, or zero signal. However, signal sources that are close to the lower coil (relative to the baseline, or separation

between coils) would couple significantly more flux into the lower coil than into the upper coil; this would result in a net flux in the gradiometer and hence the signal would be detected.




Figure 3: Response of gradient coils relative to magnetometer response (1/r<sup>3</sup> suppressed)

For objects that are close (relative to the gradiometer baseline), the gradiometer acts as a pure magnetometer, while rejecting more than 99% of the magnetic signals coming from distant objects. In essence, the gradiometer acts as a “compensated” magnetometer.

Normally, SQUID magnetometers (and gradiometers) map the axial (BZ) component of the magnetic field. Obviously, using three sensors, it is possible to monitor all three vector components of the magnetic field. Additional channels of SQUID sensors can be used to provide reference channels for electronic balancing. Portions of the reference magnetometer responses are summed electronically with the detection coil(s) output to reject common mode signals from distant noise sources. Electronic balancing can be used to create an HTS axial gradiometer from two HTS magnetometers.

Tristan and its key personnel have produced a number of measurement systems for a variety of applications. Some of them are listed here:

- Tristan has multiple single- and multi-channel SQUID magnetometers for NDE and paleoarcheology use. These are state-of-the-art systems, some with spatial resolutions approaching 1 μm.
- Multiple single- and multi-channel SQUID magnetometers for biomedical applications for animals and humans. The Ferritometer® is routinely used for clinical assessment of iron overload diseases. This system is a turnkey operation including patient scanning bed, computer control, along with complete data acquisition and analysis software.
- Tristan’s magnetometer systems are based on its iMAG® line of commercial SQUID electronics, which have been supplied worldwide to both end users and OEMs.
- Tristan’s model DRM-300 geophysical rock magnetometer uses closed cycle refrigeration to eliminate the need for liquid helium and reduce operating costs. This technology is available for use on many of Tristan’s products.

Systems built by Tristan’s present personnel during the time period of 1991-1996 include:

- A DC and AC susceptibility variable temperature and field platform. Twelve systems were made. These systems integrated SQUID magnetometers, sample motion control, sub-mK thermal control from 2 – 350 K, variable applied fields to 17 T and truly user-friendly automated control software. This product demonstrated Tristan’s ability to produce state-of-the-art complex analysis equipment with

minimal user requirements.

- A six-channel system for Vanderbilt University for general-purpose NDE studies. Comprised of a magnetometer, dewar, electronics, software and multiple magnets, this system has extremely high sensitivity (10-14 tesla) and sub-mm resolution.
- A three-channel Superconducting (SQUID) NDE system for use by a large Japanese steel company, comprising magnetometer probe, dewar, superconducting magnets, custom electronics, and custom software. Using a welding robot, this compact system is scanned over samples.
- A dual-channel magnetometer system for use by a private company to study materials for nuclear-fuel rod integrity. The package includes a magnetometer probe, dewar, computer controlled sample scanner, electronics and software.
- A compact (12") six-channel high sensitivity susceptometer capable of generating tesla fields and operating in both vertical and horizontal orientations. The ultra-compact system, when attached to the end of a robot arm, is used by a large Japanese nuclear reactor inspection company for scanning the interior of nuclear pressure vessels.
- The first commercial scanning magnetic microscope (SMM-1000) to study small electronic circuits and material samples. This comprised a dewar, cryogenic sample handling stage, magnetometer, custom software, vacuum system, and custom electronics. It is comparable to a SEM in complexity. Nine detection coils were fabricated in a linear array with 100 μm coil separation. Spatial resolution was at the μm level.
- A mixed stage (Gifford-McMahon/Joule-Thomson) cryocooler that routinely achieved 2 K.

**1. Laboratory Applications** include measurements of current, voltage, resistance, magnetization, etc. along with exotic (General Relativity, magnetic monopole) applications.

Current:	$10^{-12}$ ampere/√Hz	dc Resistance:	$10^{-12}$ Ω
Magnetic Fields:	$10^{-17}$ tesla/√Hz $10^{-14}$ volt	Mutual/Self Inductance:	$10^{-12}$ Henry
dc Voltage:		Magnetic Moment:	10-10 emu

**2. Geophysical Applications** include oil and mineral exploration, pollutant monitoring, magma flow measurements, rock magnetometry, paleo archeology, etc.



Tristan HTS SQUID gradiometer in flight

### **3. Non-Destructive Test & Evaluation (NDE)**

NDE scanning systems are used for defect detection, corrosion measurement, magnetic microscopy, etc. Some examples of SQUID NDE include:


### Intrinsic currents measurements, such as:

- Remnant magnetization
- Embedded magnetic sensors (see figure below)
- Flaw-induced perturbations
- Johnson noise in metals
- Eddy currents in an applied ac field (flaws)

### Hysteretic magnetization due to:

- cyclic stress (strain)
- simultaneous dc & ac magnetic fields

### Magnetization of paramagnetic, diamagnetic and ferromagnetic materials in dc fields.



SMM-701 magnetic field scans of an embedded strain sensor under a 4 cm thick concrete overcoating. A: Bare sensor showing dipole characteristics, B: sensor under concrete, C: bare concrete. Image D = B – C is a digital subtraction of B and C showing that it is possible to image objects deep underneath magnetically complex coverings. The scans cover a 6 cm x 6 cm area.

## 4. Medical Applications

### Studies of the Brain- Neuromagnetism

- Epilepsy
- Neonatal and prenatal Brain Disorders Presurgical Cortical Function Mapping Peripheral nerve and spinal cord studies Drug Development and Testing
- Stroke
- Alzheimer's
- Neuromuscular Disorders
- Performance Evaluation
- Animal Systems

### Studies of the Heart—Magnetocardiography

- Arrhythmia
- Heart Muscle Damage Fetal Cardiography

### Other Medical Applications

- Non-invasive in-vivo Magnetic Liver Biopsies(Ferritometry)
- Studies of the Stomach—Gastroenterology
- Intestinal and Mesenteric Ischemia

- Lung Function and Clearance Studies
- Peripheral and Single Nerve Studies
- Organ Transplant Rejection Risk
- Blood Flow Disorder



Tristan non-magnetic dewars




Model 607 biomagnetometer



iMAG® Electronics and laboratory probes



LTS SQUID sensor HTS SQUID sensor



SMM-701 NDE scanning system



Model SMM-770 Scanning SQUID Microscope



**Tristan Technologies, Inc.**

**SQUID Systems for Biomedical Research**

**Applications**

- Magnetoencephalography (MEG) for functional imaging
- Mapping GPM, events and spontaneous rhythms
- Magneto-cardiography (MCG) Adult and fetal magnetocardiography
- Spine/Peripheral Nerve Functions- measurements
- Liver Iron (hepatic iron stores) Assessment
- Intrastrial ischemia
- Magnetoencephalography (rural transport mechanisms)



INFJ 1...m inert New...

**Features**

- Safety and Convenience — noninvasive measurements
- High temporal resolution of signals
- Optimum placement, with sensor urine) to sources
- Integrated systems, ChStOre Cleweirs, aette les, patient beds,

**Technologies**

- Advanced Josephson Junction Technology
- Microprocessor-based multichannel control electronics
- Complex Custom Design Coil Arrays 2 mm to 2 CM diameter detection Pets \* magnetometers and grad ornaters (1st and 2nd order)
- Custom Dewed and Cryocotes systems
  - o adjustable tails coil-to-subject distances less Than 5 rim,
  - o ultra-low boiloff flutter; cryogen free no Gryoocolor assisted
- Custom Elect arcs and Software
  - o electronic noise cancellation hard,, and software
  - o sophisticated data equation and monitoring software
- Soptrystated Modeling
  - o finite element analysis for OM nal coil and dewer deegn
  - o thermal analysis for highest tamer performance




uu J. Sm. mum noacnr, Lworpir wompreer tah...

**Tristan Technologies, Inc.**

**General Facts about Biomagnetometer Applications**

What is measured — bi onnagnetic gelds





Etiomothoeel Applienilions

Ekinfero schtly generates son expending RImagnetic (RIM) that can be detected with Tristan SQUID magnetometers. Arise CI or MRI imaging yield struotnait infernal fiat, demeanis.

**Advantages**

SQUID based measurements are noninvasive and inner rity sale. The detection system does not contact the subject, rind subject preparation is simple and fast. Magnetic adds sutler Stub lass, OlatEken Item \*Maga potentials, pemating ORM Seelsio 14411111040n Insulating healers (skull Moon scalp, end anatomical open spaces) do not emanata or dorm magnate HAM they az with electrical potentials. Moreover, mere are magnetic mean/semens, which here n electrical analogs (e.g. stolic magnetic Melids, insomniae, susseePotilly). There era also magnetic measurements where a noninvisio procedure it 11.P.Sthht tag, total cardiography). The measurement of a magnet. Hod is a refermre-free



he 3 1010Poin Finc 'pow or Ne.NNN

- Brain mapping Magnetoencephalography ear map mean, end spodeHow nottely MOO severe nekmetets at the Volume at cedes, while distinguishing events separated by nalscercets, 111113 prostate fatigue information kcal con c&c comma rice.
- Aquil Heart Magneto-cardiography 1MCG+ uncladas dynamic metre at heart function, and can be used to nrep salroca NA rifsinglesh Werean as normal and abnormal events
- Fetal I term

The fetal MCGs a unique capability because the mos hose rafora Machi menta tion mold-no. ENG Officuf of impos-

- Sible Magnetic betas are blookod and fetal abnormalities end arriahmes can be see, as early as Del wreke el gestanon
- Spineripheral function

Naive COncarCON INVOIMS rapid axonal spikes, of a duration on the order of 1. refso Testena unique ultra Robspeed system exveceelno.oca sample, per second per channel and can detect a 2-noes change over a 5 and d&Lrnw along a seiner cord er porthema neve In


1111eStel wall necrosis is difficult to diagnose without evasive surgery. Research using glae OWI Ma OriCetes thel Scillas could replseame and esky

**Tristan Technologies, Inc.**

**SQUID Magnetometers**

**Principles of operation**

geomagnetic measurements encompass many oohs major organs alto bony, includRng heart brain. Hoer, kings end muscles Diagnostic procedures based on Etemagneasm have enlarged in cardiology end neurchey\_sad are emerging in gaStreentereley as web\_



- Detection eaglet to sense changes In enamel meg. relic 'held and Lanslorm hom roto en electrical erant
- Input mal(s) to transform Me resulting current IMO a magnetic flux in The SQUID sewer, ohotontos to transturat the applied flo into a roltage output
- A cryogenic region to keep he SQUID (normal and coils) superconducting
- Signet omsdnierirm kneurones to convent,, '...-h-he Signme ado .ollages pea n, p\*, e- normed
- Acquisition hardware end software to ours,, store and anaysis data

OS SQUIDasesteale

Pa Zol feentze s. gerff nreeprj fr, N, N, al Ntrafcre: Iga mppLNI TNINret

CV	C.. 1.04. IRni Ona. t. n l	firm, TIAlylgritnIN!
SQUID		I. caritrou
input c.ni		D a l e M j *

My 2.1. CLUv riteov FNINNow, S, Fe no UNNN JenCier

& Adam Design System sensilikty is a furionlen el SOLID serssors sersaikrity l eteollon cell diem.,

number of lures, end magnetometer/gala...1 design. Dstecton coils may be gm.

oremten rather than mairtelometers, to reject FVirdlnent none Sella. fifti, rUr of reflection channels. arangostm and spaug between chem. are all am tailored to fiser arautOROna.

Because enviermental noise vn impair det even/lot learnagnetiv signals,, systems tff

- O. acquisition and real-time menduring system
- TIM party source analysts sal... eats cornealobology with scarce analysis whams

R/ 7.0 LOAMY of enMenomeNal now

**Tristan Technologies, INC.**

**Tristan Model 637**



The Model 837 is a hay integrated magueac source repaing system with 37 sensor channel positions. Measurement coverage is nearly 500 cm<sup>2</sup>, and the array is dense ...CU to oaks imPortaint Nomadic&c maps with it single instrument placement over

\*\*\*\* or file/Tile Design: Each channel peoken wn spew no wrote from 1 to 3 (a thogonal) sensors, yielding up to 111 independent sWRces of data. The typica sensor cas &emeter is em in hesagortat dose packed among . Sens, geometries

ten be Mitered (au) changed, to follow he reesadher h adjust he spatial frequency measurement capability of he system. Systems are tailored bus, applicalaff. 1st 02Lxtemter greZernoter Optic. IE14.111 to CP. melea, and tooth of Wm...An &element reference channel system Works to recites the influence of environmental net\* SOL12.6.

System testurns: Systems typically include a pantry both e.ya. NIL end rotation) and adjustable Patent bed along with hardware and software for date acquisition and sepal analysis. The system Vatures cairnyesoman construction.

Ing 2 - Ad adyler! 3 possf, Nprom h, cone lteernt N men, 'atit' 'j, ... ehamek and N or femme, draoanely Ateddcd r i ow, and gartn, oot perm hod ars NhANNF

**Ag ralida hone**

- BMagnetoencephalography(MEG) S\*noie placement man ot sensory conical swfoe 7.1 abnormal spontaneous event Much as an Inter11 epileptic spike or skew wave)
- etwarlocardiography(MCO) areitertalMCGI YAM a gantty end real-tee data monitoring the user can adjust he system oositho ia recording has, (or fetal hewm data. The 37 channel array la suited for now swum Rhea to separate mammal from Real heart slommm.

- Ass\*ss.e. 11 of Sp.F.I.Peneh\*1qL N,, Fahation Because spinal signals involhm actin potentials, here Is en opthn to new. Man In asocess et tatopj tOrp66 per second per channel This permits daldng or ormal events, and allows tracking of vary rapid changes to oose conditioe

oww' .ki \*Ll' tINI'

f'sri nk/' niip/ wi, " wiewi

iga •ll..-oP^/ .i+V Alt.e

j

eriterAle/ Rilsteov «'11., / -)).. " er

minute, characterize intestinal signala, which are I W.2.0nm Neme, (15 rts.ter) port ITONIil radial Atm. 13 of am 37

associated with intestine! ischemia. Data acquisition r7odfun.r hope 3 ciats IVAN? vhdtruteIn.

includes nvspitu tonmorral Stan. and eme-frenency eIRays that porn& areaud mat ante mateforam



# DRM-300 Rock Magnetometer

The Trieste rock DRM-300 is a closed cycle refrigerator with a 43 um SQUID magnetometer system or measurement of magnetic susceptibility of geophysical samples. The well SQUID background noise is unparalleled sensitivity. Ginned cycle refrigerator eliminates the need to winter build helium. 16 milliliter nitrogen is needed for recovery time.

### Features:


- Three orthogonal detection coils
- SQUID detection circuitry
- 10<sup>-11</sup> Am<sup>2</sup> Sensitivity
- Wide Dynamic Range
- Room Temperature Bore
- Closed-cycle 4-Kelvin refrigeration
- Self-replenishing liquid helium ballast for Quiet Mode operation
- Automated Sample Insertion Stage
- Internal Superconducting and mu-



MM-site eon, heraweinuen (ryweitr etneetwize weid-ve)

Trieste model EIRKA 300 Rock Magnetometer offers technical enhancements to achieve superior sensitivity and dynamic range without sacrificing portability. It is used for use. Successive Quantum Interference Code. SQUID are made of doped and anisotropic magnetic moment of samples placed into the sensitive volume or the detector coil. The change in detected magnetic field is directly proportional to the net magnetic moment of the sample. The detector coils in the Yodel ORM11-300 are wound in a non-coplanar configuration to provide a range of uniform sensitivity at the center or each coil. The NM Intim separate detection coil sets are configured in a staggered fashion.

are the three orthogonal component (X, Y, Z) field



A GiliinrlfAcktah.n closed cycle refrigerator Beetles seesues return b sunray the cryogenic environment. The

dultion nuns 9x2U0 ammo. and auptoreneling Shaiti are kept at operating temperature by thermal mixing in a liquid helium ballast reservoir. The clear bore sample tube is kept at room temperature and has a narrow diameter as large as 19 mm diameter (sample size 19 mm, on either side). Ann/organic lemmenallunt oonimbrenceure mrdikaltnn subby or Ire SCIII0 sensors. The closed cycle refrigerator is mounted far from sensors to minimize vibration along the sample path.

For the ultimate sensitivity the DRM-300 can operate with the cryocooler tied off for periods of more than 10 days. An optional vibration isolation system is available for use in continuous operation.

Further customization and enhancement for the DRM-300 is possible through the offered options.


### Options

Further customization and enhancement for the DRM-300 is possible through the offered options.

- **Vibration isolation system**  
The DRM-300 is designed to be used in the center of a building. The vibration isolation system is kept at operating temperature by thermal contact to the liquid helium reservoir. Additional vibration isolation system is available for use in continuous operation. The independent detection coils and weighted magnetic moment of the detector coils are mounted on independent vibration isolation supports which rest directly on the floor, independent of the rest of the system.
- **External Magnetic Shield options**  
An optional mu-metal shield mounted outside the bore is catered for further reduction of external magnetic fields. Frisian man also supply magnetically shielded rooms or limeaxis cancellation.

### System Specifications

<b>SENSOR:</b>	Three Superconducting QUANTUM
	Interlarsense DOViceS (Slaters) operating at 4 K
<b>SENSITIVITY:</b>	10 <sup>-11</sup> Am <sup>2</sup> /NH, ennu white anise
<b>DYNAMIC RANGE:</b>	10 <sup>1</sup> Am <sup>2</sup> (197 dM, higher ranges available)
<b>CRYOGEN-RI-M5C</b>	Not needed - Self-replenishing liquid helium ballast for Quiet Mode operation
<b>HOLD TIME</b>	Infinite, 2r days Mb Cry/Wader oft



# Model SMM-401 nanoSQUID

scanning microscopy for imaging magnetic nanostructures

- SMM-401 uses a superconducting (50,14) sensor to provide outstanding spatial resolution and high sensitivity

### Features:

- 100 pm spatial resolution
- 1.4 pT/Hz field sensitivity
- Room temperature sample
- 25 um gap between sensor and sample
- Non-magnetic scanning stage
- Low helium consumption

### Applications

The SMM-401 is particularly desirable in the areas where high sensitivities. Especially at low frequencies, a magnetic field of 10 pT is required for the detection of magnetic nanostructures.

When imaging a homogeneous magnetic field, the SMM-401 provides a resolution of 100 pm.


Resolution: 0.120 Goan of Pointentesor

ter, ei

Magnetic microscopy image of a 2 um by 1 um slice

of a martian meteorite ALM400Y overlaid on top of a dark-field photo of the same slice. The green shows the field intensity.


The image shows the magnetic field intensity overlaid on a dark-field photo of the same slice. The green shows the field intensity.



0.120 Goan of Pointentesor

ter, ei


Magnetic microscopy image of a 2 um by 1 um slice



# SQUID Magnetic Scanner

The Tristan model SMM-601 Magnetic Scanner is designed to measure magnetic fields with a spatial resolution better than 300  $\mu\text{m}$ . It can be used to image diverse objects such as:

- subsurface cracks and flaws
- embedded magnetic sensors
- composite structures
- corrosion sites – hidden or exposed
- impurities in metals and



Scans of 3.5 and 10 mm holes in a metal wafer


Besides measuring magnetic fields, the SMM-601 can also be configured to detect:

- + induced magnetization
- 4. aging and stress in ferromagnetic materials
- magnetic susceptibility
- eddy currents
- magnetic hysteresis
- Barkhausen effect




Mode sAverli Scanning Sales In microscope with dc endure big na magnets

- Tiro SM151-501 in a fatly featured Measurement system that ,Noise the merle mama inagenc image of the Abject being measured over the untine In — 50 kHz frequency range The adjustable rail viewer the spacing between the detection coil end sample N.) to be his small as 3 turn.
- Ice lima frequency Iresponae means large penetration depths (deep penetration). Another advantage of the model SWS1560 i is its ability 10 Oidiele hi rEin fields. This alleles il 10 Matt SOSCLMignity measurem.rs on the same sub.mmm spatial resolution scale. In addition, it can operate in ac with do—en td 1n bandwidths for eddy current measurements.
- The shirf eel allows computer controlled scans of objects over a lmbhavi (1 K 15 mu) ores with 25 pm slapping amiability
- The Use of a dc SQUIT) sensor gives it opatelleled sensitivity Its cal phase response allows bola in phase and quadmrture information to be obtained Swithout distention. If ultimate sensitivity is needed, larger demotion rolls with resolutions exceeding 5 frin/iz are available. Additional detection coils can be supplied to give vector inbrma lon The Shil-601 requires mdmal seep. Automated setup and computer



glannae add moans of an oneadand strain serisee under at em thick caner, oVerco, A - barensene, showing uharadvinathys- sensor rider cow., C - berm umcrate. Imsgo 0 = 0 — Cis dieflal subtraction of D and C chovin taint it la coala to lmas olibita deen enkn 114



The stendeved modal eht'ol-801 is ,margured to detect electric cotton, and to measure morlan magnetic Odds. It includes a Single.Chunnei scanning SQUID Maamlomemtar Probe. INVAie 551110 Electronics Cristogenie dew, Room Temperature Scanning Stage co-make contiol ore Mtn 011pdeller System, end imaging Sothic. The model Stall-edl ran ha sumvimen Mill adeoneel capawites th extend its measurement capabilities.

### OPTIONS AND ACCESSORIES

Aditidonal DalooNan Chomais: The model SHM-Ell's measurement capabilities CIENI be extended le multi-chennet oapabilates. This can mean either vector IS, & end By) wpm/Wiles .ar

ndMionW Verli.1 .....re wt. to reduce rmasomMent Coins lodtallan Ch. Meg even NSO tie Odd. fa 60EG Wl,m anterantennal noise In OkapISSIVIL

Wadicahae rotnrl.almaap linrlfmlnc

Seen Aram Lerner scan areas and higher mann, stepping 125 pni sternard) are available per mowed,

ca FMS Capabily. This option mows a small an t'gellnefic adde le exceeded el same. Tire Fell is vertical ,su and can have a peak-teapeah magnitude at La 1 oersted This cerrorily is Cr Whaler interest eyes eddy careni measurements eve desired

NW/A1W ALPS(B, and 00 Sham Ingham: A lonvzontial field sheet inducer mach ran apple an ec magnetic field parallel in the teat surface In in.ce a large intended eddy currant in e desired onteanent. can be used to image cracks cr material Ion Jeep In eunduebee tug., aulmbnbnr sthichnuies.

as Field Componsallon Ebotranica. Men en ac mon, is dimity coupled into the syglern, the result. signal learn, field colts) may be trot it lamer than the atonal hom the sample. In the case of a logourenette materials sunk so carbon steel, The induced MarrellValieg (Man Wen n 110011 01 WO Matt in pulp lame and me ripe., range or he dextraumen, mishear Thar not he elmaysto le track this large signal ahle 101 Tasokano the small Nona Awn def. mine metal.

To rininnoe thus, Tu. n. t. tub4y an Jr Compete ellen system to nil Pie as signal In the detection era and eren the igrluctual Venal rit the chest hang me.ured

### SPECIFICATIONS

srAisoR, Lien teissarhehlo suraoamossina woken leustrareaso d.v.ta lScesto

SPATIAL RESOLUTION: Getter th. 300 pin

SENSITIVITY, a 10<sup>-10</sup> haelehgte /50 ITivhgl rov 3 mince.

DISTANCE TO SAMPLE Adinet&h lo 0.1.5<sup>mm</sup> A MM

OPERATING BANEAOMOTH: da - 50 kHz. Measmerrients net be

### TRISTAN TECHNOLOGIES

61 IA COMehelicric 00Crl 001, 0010  
105 Say, dug, ca 52121  
55810511-2706(thg550-2799

# Scanning Magnetic Microscope

The Tristan model SMM-770 Scanning Magnetic Microscope is designed to measure magnetic fields above a planar surface with unparalleled spatial resolutions,

Using a liquid nitrogen SQUIB sensor, it can be used to image room temperature objects such as:

- + traces on a circuit board or multi-chip module
- shorts to ground planes



Ma01..11nagc Ordollar

- insulators, ferrous and non-ferrous metals to detect cracks, voids and corrosion
- nanoparticle distributions
- flux-motion in HTS materials


The SMM-770 can also be configured to detect:

- induced magnetization
- magnetic susceptibility
- eddy currents
- magnetic hysteresis
- micropaleontology



SUMO Microscope Loll o: a nd ac have mans.

- The Shilv-770 is a fully lectured measurement system hal allows the user 10 extract a magnetic image of the object being measured overlie entire do —10 kHz Erequency range.
- Whether the samples are circuit boards, multichip modules, steel or aluminum plates, composites or even plastics, the SUM 770 cen measure surface and even deeply embedded sources with a spatial resolution down te 50 strn.
- The use of a High Temperature Superconducting do SQUID sensor gives it unparalleled sensitivity with the ability to measure fields smaller than 20 pT/Uftz. Tristan's FITS sensors can also operate in applied magnetic fields up to 10W oersted.
- The Shl M-770 allows computer controlled scans of objects over a large (15 a 15 ch/n area frith 25 Pat stepping capability with sub4incron stepping available-1110user Fee the ability to preprogram the wen ceardinetes.
- The SMM-770 (emigres minimal setup. Automated setup and computer couglal makes Theasurements rapid and repeat able. System sort-were provides the ability to control the critical system components, acquire data from be SOL(111) lessor( and analyze illa data no determine the magnetic progenies oitile sample being ineasured. The use of open architecture software allows the user to cusinthaee nearly all



an-juo. of Ma Medal S.11.77 system Naucbaq, n1 II accessories

The Areidand Shh4770 is configured to detect elecui entrants and 1111.5011 remnant magnatik Voids. It includes a Shih-Channel Scanning Shim telegenometer Preto. MAC. 501110 Iactronics, Cryogenic Vegan. Roan Temporal. Lminning Stage. Computer Chantal and Cate Acquisitive System, and Imaging SIAVAIR. The

Stas.770. can be supplied Yeah addittatit capabilities b Whigil 110 1710.nriarriant haepOnlatlei

### OPTIONS AND ACCESSORIES

SCAN AREA: Latger.en are. 30 cm a 30 cm) and tighthe lesslbnl atimp l. (20 um stranded) no ovulaf upon mallisel.

Sularlrinten OF 1051-1 SQUID Salmon For me.unnals Maim .1sittvitas significantly Pak. 20 pThilde are needed and ultimate spatial naselliton is nal SO Important. Atkin can substitute the madgthHrhals.or with S sijn&aney law rs semi. detecIsie "al

with a .thsPivira d.c.r than 10 olio-lo 11TMA1<sup>3</sup> mm carmen with a sensitivy holler than One paltht teal:30.711M.

Woricanuhartheructelunimaiso mat

Norinmtal Flak, DA add 13a Mae, increco, A horizontal Sail sleet inducer. With can apply an . marmot field parallel to the Wt mitre b induct s tinge extended eddy current in a remind enentailin. can he wed to inage clacks or malc..1 b. deep in conductive Mg, aluruntirell structures.

as Rent Compensation Eelektrnkr When imaging, conductive Materials. il on an signal is dinctly coupled into the system the memo signal tune the held 101151 mss he most lairs then lee signal lam the sensA. To minimize this, Tristan can supply an ac Co mprearg n system In nauthe no strollin the dotageen gni and

etrant the \*0150 mgnal In the object being trairstained

### SPECIFICATIONS


SENSOR 1-1gh temperature superconducting quantum arlorerenco diets 1500101 op-ember; at 77 K

SPATIALRESOLUTIONVReiter than 51 um

SENSMATTV, 2th 10' trantIA He (20 Tsdsstri

OPERATING SANDAIDT/M: dc - 10 kHz. Measmerments can lie made al any inagya nay, ganderatris show 20 MFR am available.

CRYOGEN/O COOLING: To ahead bw tregueirey vett balms 20D hg the system uses IWul narogan b cool. sen.r



Scans of 1, 3.5, and 10 mm holes in a steel plate

### TRISTAN

III TECHNOLOGIES et.,  
Cwrea brew Court East, Salle  
106 San Ciao, GA 92121



## Magnetometers for Geophysics

Mineral surveys, magnetotellurics, magnetic detection of induced polarization, and other magnetic detection methods are important geophysical tools. Superconducting magnetometers and gradiometers offer several advantages over other detectors commonly used for such measurements.

- Constant Sensitivity from dc to 10 kHz
- Magnetic Field Resolution of  $10^{-14}$  Tesla
- Gradient Resolution of  $10^{-1}$  Tesla/meter
- True dc Response
- Flat Phase Response
- Wide Dynamic Range

Tristan manufactures the most complete line of ultrasensitive geomagnetic measurement systems available. From compact single and three channel Magnetometers to 8-channel tensor arrays, Tristan offers a variety of fully configured system packages for geophysical measurements.

The basic geophysical measurement system offered by Tristan is the model G377. It measures all three vector components of the Earth's magnetic field ( $B_x$ ,  $B_y$ ,  $B_z$ ). The small size and portability of the model G377 makes it convenient for field use. It can also be supplied with different size dewars for airborne (model NLD-530 dewar and borehole (Model NGD-530 debar) use Planar Gradiometers can also be



Tristan offers three basic sensors for geophysical measurements, the FITM-5 and the higher sensitivity HIM 16 are magnetometers (13, Br, Be)', the optional I-ITG.10B measures planar gradients (d13, Ids).

The model G377 can be supplied with fewer sensors or a mixture of magnetometers and planar gradiometers if needed. The picture below shows a single channel planar radiometer (I-ITG-10R sensor) being used in airborne measurements.



Model 701G system mounted at end of boom

Model G377

Model G377

Operation Principle: 3-Axis 77 kelvin do SCIL111:: Magnetometer - Measuring the relative change in magnetic field simultaneously in O, B, and 13, BOORS.

Range: 7'5 pTHz

serheavkffh: do Iota kHz

wider bandwidths available

Slow Rate: >1 pTrsec speak to peek)

Sensitivity: 50 TIT7VHz, HTM-13

20 frNW, HTM-16

1 frmHz: HTG-10R

Cryogen: Liquid Nitrogen

Volume: 7 liters

Hold time: 2-3 weeks

Power: 120 or 240 VAO, 50 Watts

112 Volt Battery Supply Optional)

Outputs: Analog, RS-232 or IEEE-4138 Visual Alphanumeric display Controller. 321 MM 0 121 nrit X 300 mm (12.6' wide, 4.8" high, 11.8" deep). Weight: 3.5 kg (8 lbs )

NG13-1030 dewar (7 Hers) Standard on G377 406 rem high, 250 mm diameter OR high, 9.8" diameter)

## T877 Tensor Gradiometer

TRISTAN TECHNOLOGIES



The Tristan model T877 SQUID tensor gradiometer is designed to measure magnetic fields and gradients for geophysical measurements.

It is a valuable tool for:

- Magnetotellurics
- Controlled Source Measurements
- Borehole Measurements
- Transient Electromagnetic Measurements (TEM)
- Unexploded Ordinance (LAO)
- 4 Magnetic Anomaly Detection
- + Environmental Waste Detection

Superconducting magnetometers and gradiometers offer several advantages. Over other detectors commonly used for Magnetic Anomaly Detection, MagnetoTellurics, magnetic induction of induced polarization, and other geophysical 11.15mEntenle SLVeronntrmln detectors' offer constant sensitivity from 1n to 1000 of kHz (or higher),

tad magnetic field resolution op nT41z with magnetic' gtedleal resolution up to le ISPT4HA and a dynamic range of 40 dB. Dame systems are well 11.0100 10 illea oe, be1ne ilahweight, [suabe, rag to set up, and easy to or.

The T877 masnetometerradiometr offers owsl important advantages over other magnetometers. It is a vector magnetometer, in contrast to the proton mensation dem= whwh responds only in the magnitude of the Seat With a three-axis vector magnetometer, both the magnitude and direction of the field can be determined. With eight sensing elements in s tensor configuration, the complete magnetic field gradient can be determined. Its performance is not impaired by the presence of large Anieries and — unlike Niagara devices — SQUID magnetometers do not saner.. In comparison to large induction the itor in col awkward or cumbersome Sr deployment and use The T877's response avoids giving andue emphasis to high frequency phenomena such as the ubiquitous lightning induced stenos.

phaseahfts.

The Tristan Model T877 is a tiott-prmo0 reglad, highly sensitive superconductin SQUID maganmcieigradiamter design: for geophysical exploration and measurement. With the full tensor configuration, it iv possible to doth complete clunacter.nan of may' erie dipole sources Sr long nags obtaining kmahmation and clualitication information. MIE has been shown theoretically by Wynn and demonstrated in rle Field All pat necessary is knowledge of the magnetic field compe.. (Hr, Hr, Hy) and m< live unique field gradients (a3agg, PH, 04',

70s0S, Mies) The Titt combines etgkl individual magnetometers into anarray that yields all mummy Odd endTarr BMW, twang gradient components.

For urhosne operation, Tristan can supply custom-dewars-including-herisontat or other Wynn, a se', Adamed superandganmeramotombantansrAnys mats horyosgat praegsprnTectis.s. rgrTrans out MasinOt, 071-7711M31



### SENSOR CONFIGURATION

• The /neg... field vector,  $\vec{f}$ , can be expressed he cones of Cartesian CONFOR15 //-. For each component, there are three sparial derivatives along orthogonal directions, generating nine components of the non rank magnetic veld gradient Sera This two can be represented by the Finns.


$$\begin{aligned} &OH, 614. 58 \quad ag, \\ & \rightarrow 7 c ; \\ & idr OIL \quad ahc OR, \\ & ) \quad 7.7 F., \end{aligned}$$

• According to Maxwell's equations, only five of these tensor elements are independent which is the case for the SQUID tensor gradiometer.

$$\begin{aligned} &IX \bullet 1.11- \\ & 0 2 \\ & .H=Z(2+2'I \end{aligned}$$

oar of tikm hm,11 we

d e a d ,



Data from Controlled Source Measurement cal a 10 n diameter hollow MO me Mica copper sphere showing data horn the foe indoPendont tensor gradients Dale was collectory ea ow Intervals.

### SPECIFICATIONS

sENS'On', High temperature superconducting quantum intarlenrsm Cna (SCU)Etst operates 017711 OPERATING RANGE: t 900 nT

DEWAR: neters  
VOLUME: HOLD  
TIME POWER, nominally 2 weeks

129 07240 V... 50 Watts Ito Volt Battery Supply Optional)

CONTROLLER: Analog 0.3 Volts

- Model H7148HTS or SQUID Magnmemmer Sensors (S) DEWAR
- Medial NGD.1099 Liquid Nitrogen Dewar avowereu dm 1,0131316 dams e.11.able on 11.068
- Model NO1-100 cryogenic insert and oriet:genre robins
- Model MC-393 IMATS01.11 3 Electronics COMM Soil
- Model iFk-301-iIFkre-Locked impis (r1)
- Model CC-09 six meter lIberopto composes cables (S)
- Masal and accessory perk lira. on

55232 or IEEE-436 41.1.1Alpha,rtleato dienNy 521 rrm wide, 121 ram Nutt KO rum deep 1125fx 4.6'0 11 8 0 9.13 kg (5155)

TRISTAN TECHNOLOGIES

ihSSComerstonn Court Ent, 5,01105 San 0: 0e. G. 92121

## Biomagnetic Liver Susceptometer



## 0.142 Biopsy

The Biomagnetic Liver Susceptometer is a diagnostic instrument which measures iron stores rapidly and non-invasively. Its advanced design with a superconducting magnet and SQUID detection system gives an accurate measurement of iron concentration in the liver and spleen for adults and children.

- Non-invasive  
Replaces Surgical Biopsy, for Iron Measurements
- Eliminates Discomfort and Risk
- Allows Pediatric Measurements
- Direct Measurement Method  
Accurate and Reproducible

## BLS Methodology

1.01' Snsseptometry (BLS) ekulorts the effects. of magnousin and supacone-heus in Biological ma ten Ms such as rennin and honosidenn as mallsartrauled to an molled magnetic field Maracsagedic behinor' while untr and hod. tissue are,em weakl, topched (thannmehel Ferromagnetic Platenals c g, nickel and anal. are shong/4 nthacted n, polled fields No Bengali ty penurring blane Bunn la fmn:triansatin

In Is BLS method. a weak magmas, NM] n#0 20 Ind litesla is 'enoned tathin the bed), borate, an external supewconduchnn field manner similar In thra oscd ins MRI scone, bat a hundred trees weaker.The applied fields me =elated M a superconducting magnetometer known as a SQUID (SimeremithaturnQrantsint Interference Der lcc) The <sup>59</sup>Fe D. raing, stom his ill. ability tip F.ASFFIFFFFIMM W' niamicnt nold al the pax no billion heel.

Wwi, an organ, such as the hoer IF placed ina mannenc field. Ir mall slight', dsarcot the aughed field, hole (Fuer a normal' or samara. the kraa hetet 1611 no redoned slightly. If drainer is lmo

orecloudad, the local field will be enhanced. Hem the pliangt In **Specifications** Magnetic Yield 20 InT el cod Dec. stab& n better than u.l pphndour. 5 Gauss bins Il am rum Jut., lad Detection Coil Sensitivity Int 174117

Liquid Helium Connety 35 1 lers11 ft l., hpid ham

Parkin Bed Capacity 131 Lg.

xia min resolution

Daze Base. Open file stmctbre with abdpd Mauston™. to user pref, reases

\*Total System Noise. 00.02 meg \*Fel coucknknwOor (met w111111. mansinvd ia 'her phantom el a distance OF 15 rum llsonment Range 115-30 niwg lwer ucighfl, ari, Diem 10 to <sup>200</sup> mgia (do weight,

## Options

For researchers dmerasted to dvanding measurement capotaline s\_ Trsten offers the following options-

- Addition, " Detecion channel fine Inekudm a third delecion EMI with dffinent sarrisl 5S11011,1
- Atti. NOPISF C-annallition Needed For sites with high ent, minnernal poise
- All Fidel Capabilin tie-weld modulance of the mannehe field. This can ado. .t.search on alternate malhods of BLS

11. lryaauvEnx.reronrL LineSer,inner, etpalfiedaracionwaigariaeolr,

mmr,m,r-ronLienr.r.m'4, easoord.,

and, a, sme, n,1\*\* \*11., lo0ol  
, Nowee.nre n,whiti n, ab, welta, el wee e

On detected magnie Enid is 501,0010 mlakd to the troll a:neon moon in the liver

To mina nice the body's contribution In 6c-distortion In magnetic held, a eroall bag of warm is placed hizr:oh the delectaur nod dun susfrat. Since the smocedhdillyofbad, 01000 - close to Mal of halm', the resultant aoommremet a msen alt. that of a magnetized 6,1 for spleen/ moving in a uninnenn field Within i

lthamattle,jatrarvrirpnmnt, 410 -only chump Span by the &nation ends is due to no: liver lot shipentl itself

For higher accuram. our software removes the annul eon ribunon of nralyng tissues (skin bone, -MC ic, iha, Mc). This gives the line concentration at the lnes for spleen) alone, allowing ac,mrthe Moasaran lints for obese patients and normal patients 00,11011 1,1011 lir mhpken depths

Tr, sads die BLS malnV11am rem: applied up orzans. 5red i

Inas and enlamot spleens f+300 nith a total orm of

IFel = 0115.11 mein tissue (wet weight) Repeatabi ison al measurements o, cr Max ,hnds/ on stogle sbibxets -of betor Mian

OS <sup>50</sup>, has b, undeinonstru.d.



## Site Requirements

The model BLS reamims & Minimum 3.1mvr M1 nix in

r ITx no Torn) or amr n eight ee 1. 5110 kg Pkure tegaltment is 7 hVA, , tibration free 111110001 for the gam) or required and the system should he sued rp a magnacnally

quiet cmiwormoni Cont. Tn elan fo discuss site servers for magneto and labwatton mosque. Ms. Al Tristan', Enritrais me colcoed kn. a 1 ear nazranj. Some contracts may be 511F, heard a pm ide past nurruntL, esmage

**Clinical Relevance**

The standard measurement of iron stores has moved a step or needle in accuracy. This method requires a physician's analysis. With its associated time delay in obtaining results, it also assumes that trouts are evenly distributed throughout the body. In addition, the needle biopsy procedure is uncomfortable and, in some cases, significant. The

serum ferritin measurement, which has been shown to be a poor predictor of actual iron stores, is non-correlative. In addition, the use of <sup>59</sup>Fe-erythrocyte labeling in patients with iron deficiency anemia (IDA) is not recommended. The use of <sup>59</sup>Fe-erythrocyte labeling in patients with IDA is not recommended. The use of <sup>59</sup>Fe-erythrocyte labeling in patients with IDA is not recommended.



**Measurement Protocol**

The first step is to ensure the patient is in a supine position. The patient should be lying flat on their back. The depth and shape of the titer No spleen Purified

measured by ultrasound and entered into the patient data base

The patient is then positioned on a specialized couch




The bed is elevated until the pattern of the detector and the Haler hell." is tilted

The patient is then automatically lowered about 10 seconds. As the hill 101 APTS water bows 1010111e hdows keeping We space Seto eon the pal., and detector tilted

The detector is then moved to the field measured. The distance of the liner from the detector is coded

The computer immediately analyzes the data and gives a result as soon as the bed stops

The measurement is generally repeated 3 or more times to improve accuracy



With the shift 1.0 Use into account the contribution of overlying tissues. BLS measurements can be used in adults and children who have wide

variations in organ depth and CM Another advantage of BLS is that it is a non-invasive technique, giving an average iron concentration measured over many milliliters of organ tissue. It is accurate, portable, and safe. A small needle is used to measure iron stores. A small amount of tissue can be removed, and the amount of iron stores can be measured.

**Applications**


The most common application of BLS is in the diagnosis of iron overload. BLS is also used in the diagnosis of iron deficiency. BLS is also used in the diagnosis of iron deficiency.

- Monitoring iron overload in patients with transfusion-dependent anemia (TDA) and sickle cell disease, or other transfusion-dependent anemia, for diagnosis or assessment of response to therapy and iron chelation
- Assessment of iron overload in patients with hemochromatosis and assessment of the degree of iron overload on iron chelation, hemochromatosis

- Assessing iron overload in patients with hemochromatosis and assessment of the degree of iron overload on iron chelation, hemochromatosis
- Assessment of iron overload in patients with hemochromatosis and assessment of the degree of iron overload on iron chelation, hemochromatosis

study

- Diagnosis and assessment of the degree of iron overload on iron chelation, hemochromatosis



## SYSTEMS

Tristan offers a variety

of fully configured system packages based on the IMIAG series or SQUID components. These range from basic single Fenannat magnetometer systems to instruments for specific applications. They include systems for biomagnesern, geophysical exploration, nondestructive testing of materials, magnetic microscopy and studies of rock magnetism. For applicettinne that require applied fields, Tristan can supply persistent superconducting magna e, permanent magnet structures with custom-designed field profile shapes and built-in copper magnets for ac fields. Tristan's SQUIDS are available in both high temperature. IIHTS1 77 K and low temperature. [ILTS] 4.2K veratione. Standard product data sheets and application sheets are available to merry mr /mos. complete systems. Contact your Tristan products representative with your Defactia

a Laboratory Applications a  
BiOmagnatio Measurements a  
Geophysical Exploration te  
Non-Destructive Evaluation

Y Magnetic Microscopy e  
Custom SQUID Systems

The basic SQUID system consists of an input circuit connected to a SQUID




Fig. 1. Typical magnetometer system

Tristan offers complete systems or individual components, according to your needs. Tristan also supplies the Isaac components that can be combined to form the teas of a SQUID measurement system. Specific links-malids on individual components can be found on their respective data sheets.

### SQUIDS

- Model LS0120 LTS do SQUID Sensor
- Model FITM.100 HTS Magnetometer
- Model 1'170-100 HTS Gradarreler
- Model H7.100 HTS miriMAG

### PROBES

- Model SP Standard Cryogenic Cable
- Model RMP External Feedback Probe
- Model MFP Multi FdraCtion Probe
- NLI series of dewier inserts for FITS SQUID sensors

ELECTRONICS

## TRISTAN LABORATORY SYSTEMS

Tristan offers Tie 0051 complete line of SQUID mesura- ment systems available. These systems can be combined with either user- or Trietan-supplied cryogenics to give you the most versatile measurement capabilities possible.

For laboratory applications, the LTS SQUID system can be configured to measure a wide variety of electromagnese signals. VETS SQUIDS are available as pure magnetometers and planar gradiometers. Typical sensitivities that can be achieved with Tristan SQUID systems are listed below:

- a) Current: 10<sup>-6</sup> re amplitds
- b) Magnetic Fields: 10<sup>-6</sup> re lesiWHz
- c) do Voltage: 10<sup>-6</sup> re volt
- d) do Resistance: 10<sup>-6</sup> to i1
- e) inductance: 1(-4 henry
- l) Magnetic Moment: 10<sup>-6</sup> emu

Model ISMS Haan Measuring Systems: The Model eMS-H is a HIS SQUID system capable of measuring magnetic Seas approaching 30 lemtotesiarefle (1 ft -4 10<sup>-4</sup> testa). Typically, this system is used in conjunction with a 411.0 series Devar.The BINS-H can also be supplied with a planar grad. above

The Model SMS-L s a LTS SQUID system capable of measuring small electric currents with a better than War 7 X 10 re ampereNFla. With a simple pickup mil, it also can be used far the detection of magnetic fields as small as 1 fT.

Model PMS Picovolt Measuring System: This cryogenic dc velleage amplifier with a gan or 10<sup>4</sup> and a rms noise of less than 10<sup>-6</sup> volts4H2 is used for measurements of very small voltages and resistances.

Model MPS Multi-Purpose Measurement System: This system is a low impedance ac bridge system for extremely sensitive resistance and inductance measurements. Resolutions of 10<sup>-10</sup> in ohm and 10<sup>-10</sup>ts henry are readily obtained. The Model MPS also has the combined capabilities of

## TRISTAN

### MAGNETOMETERS

For measurements of external magnetic fields, Tristan offers both liquid helium and liquid nitrogen SQUID measurement systems. Senes 500 LTS systems are designed for the researcher sac desires urinate performance from a low to med. um Manna yours SQUID magnetometer or gradameter system. The series 700 HIS magnetometers offer researchers

interested in HTS (nitrogen) SQUIDS a number of convenient platforms lo pemrn magnetic measurements.

node tyre eereet orienislon

808	170	Out	1.1, 7, 14, 0.19, 8,	dst)
12	LTS	inn		
01				<90 W/n stallfarinst
cos	SITS	3	Pstif 10'11	e1 RNHo

Web the use of discrete detection circuits. Tristan LTS SQUID systems can operate in magnetic holds exceeding 9 rests and sample temperatures ranging from crIk to wall Devar. The BINS-H can also be supplied with a planar grad. above MOW temperature. Tristan HTS SQUIDS can operate in fields that can exceed 0.1 testa

### TRISTAN CUSTOM SQUID SYSTEMS

Tristan has supplied a wide some/ of unique SQUID-based instrumentation for Laboratory, Biomagnetic, Geophysical, and Non-Destructive Evaluation (NDE) measurements II your needs are unique, canted us to discuss year particular requirements. Tristan's scientists and engineers 20+ years of experience and an ever-increasing quest for refinement of its products, secures that Tristan can manufacture the ideal SQUID system to suit your needs.

SOSCR.15-nat. caress Azrue.

riff) TECHNOLOGIES

## CRYOGENIC PROBES for the laboratory

### FEATURES

- 4 Easy to install
- Multiple Measurement Capability
- Irregularian cr Vacuum Operation

### Tristan

manufactures three nams Scil J10 probe. for general

Such

- Magnetic Fields end Raid OraLliente u Static Magnetic Moment and Susceptibility Electric and Magnetic Fluctuations


Tr Intents cryogenic probes and cables are the heart of any SQUID hosed measurement system. They provide a flexible Fansmission line funning tram room temperature to either 4 K or 77 K with plug-in connectors at each end. Without restrictions of a rigid probe, a variety of installation options are available.

In all Tristan probes, construction matertals are non-magnetic and carefully seleded to minimize conduction or heal into the cryogenic Path. All probes are shielded against d intorlerence and other etectiel transients het may anent the SQUID operation. A mom temperature 0-ring seal al. lows Rowed decor operation Probes are available separately for up. grading older SQUID systems or for expanding the capabilities of a more r. AnAtky omhshsed system.

The Model SP Cryogenic Canis is the probe of choice ler simple sumori and megnetc: telliMeasremen16- Seed with MDModel 05009 low lem-

Perdu re ILTS)de sensor, measurements shown In Fig la & tb are possible. Used with the Model HTM+100 high temperature (FITS) do SQUID sensor, measurement configurations shown in rig it are possible.

The Model MFP Multi-Fraction Probe is the most versatile LTS SQUID vacuum sealed housing. The Imper LSQ201ASQUID sensor housing accommodates a cryogenic nominal hoard Oath the testleze. And nu dual mdtance Stan dards required for voltage and impedano measure. nerds. Its great versatility makes this the sonetend. ed probe when a variety of applicaloes eneequired.



the Modal RMP is dcaligred ell am moaner...nests (Fig. le \*5 and cantorelons ra nrdnks external feedback


et MO

11)

leMellif art., I WW1 tr<sup>21</sup>

### TECHNICAL FEATURES

- Model SP: Waiving Tenwensture) 0 - 77 X ISnows dependent)
- Model AMP; Melting Temperature: 0 - 7 K ELSI7IOVM weeps only)
- Standard Mutual knuCntandit 0.0 uH (nominal)
- Input impedance capacallva al nontaro frecuencies wall Z is 1120 / to
- Current Leads: d decousted fraung palm masimum Current 0.5 Amperes
- Model MS.: Working Temperature: 0 - 7 K (LSO/OrM sensor only)



III) TRISTAN TECHNOLOGIES

6105 Comerone Caw East Simle 106 Sim Dego, C, 1 93,21

page AM. 559-700 Etr. (056107 Z 90) #1144444444



# SQUID iMAG

## SQUID ELECTRONICS

- Peatneo:
- S Easy Setup
- IP Manual end
- Auto Tuning of All SQUID Parameters

### 1- Multichannel Capabilities

### 11. Single Controller for LTS and HTS SQUIDS

Tristan's MAAS SQUID electronics have been designed for the user who wants performance and flexibility. edicroprocessor-ekden hierarchical front panel menus allow fast setup for both LTS and HTS SQUID mesure- s. Multiple stew rates, gains and bandwidths allow the user to fine one the measurement process. Individual tuning of each channel ovids optimum performance in multichannel configurations. When you need the best in SQUID electronics, look ter thp iMAG series to satisfy your needs.

Thou Model INIC-373 IMAC SQUID controller forms the basis of a powerful and flexible measurement system. Its three channel capability accommodates nearly at laboratory SQUID applications without incurring the cost or complexity of eight-channel designs. A unique feature of the Tristan controller is its ability to simultaneously control both LTS and HTS devices. For the experienced user, the Tristan Multichannel Controller of-feria complete manual control of all SQUID parameters, including ties led el, modulation amplitude, "skew" level, de flux level in the SQUID (offset), heater and integrator reset. All parameters are easily adjusted using the membrane keypad and a convenient menu-driven intedece Liners who want a fully automated system will use the one-touch tuning capability that rapidly and reliably optimize the level at all critical parameters.

High-resolution ND convertors and the standard IEEE-428 bus make the IMAG controller ideal for use with computerized data acquisition. Use the rear-panel 050 connectors to monitor the high-level analog outputs. A "fourth channel" input allows you to synchronously digitize your own Mg-nal along with the three SQUID signals using the controllers internal AA) converter. LabViewtu software drivers are also available-

### FLUX-LOCKED LOOP

MAD FLLs are offered in both FITS and LTS versions. The LTD version uses an advanced bias reversal technique [hal effectively reduces low-frequency noise in HTS SQUIDS without introducing reties spikes in the output spectrum. The less-expansive <sup>LTS</sup> FLLs provide slightly higher ire-quarry response.

The Model IFL-301 series IMO flux-locked loops (FLLs) provide superior performance under a wide range of operating conditions. The Tristan design locates the FLL as close as practical to the SQUID sensors and eliminates the need to run low-level or high-frequency loads over long distances. A shad Cable connects the FLL to the probe or cryogenic cable, allowing the compact FLL to be conveniently mounted near the dewar, hut net of the way of the liquid cryogen transfers Connection to the WIC-303 Controller is via a composite cable.

### COMPOSITE CABLE

Thetane advanced design provides superior radiefrequency (M) rejection and allows for long cable runs, even in hostile environments. It is a simple matter to locate the FLL inside a shielded room and operate it using an IMAG Multichannel Controller located outside the room.

The connection between the controller and flux-locked loop(s) M via the CC Series composite cables. Low level dc power one the high-level ens-log

## IMC-303 SQUID CONTROLLER



### INIAS Controller Nadel 00.3431

Number of Channels: 3 SQUID channels Net interface to both HTS and LTS Flux Looked-Loops (FLLs)Thecontmler can operate any combination of LTS or HIS SQUIDS simul. laneously using the appropriate FLLs. An auxiliary channel is stippled for synchronous data acquisition (see below)

Auxiliary I/O: One auxiliary analog input (10 Oa impedance 50 kHz 0W) is prodded for 16-91 digitized 01 a user supplied signal for synchronous acquisition or event trigger ing. Maximum output signal is 4.5 V FS. Gain is selectable to be x1, X2, X5, 7010, 020 or X50.

User Interface: Interactive user Interlace vie large LCD display end membrane keypad. Special fudion keys and menu-driven software provide friendly operating and setup environment.

Remote Interlaces: Both IEEE-4B9 and RS-232 remote control interlaces am standard. All control sethings may be input and all instrument data may be output via these interfaces. Total maximum data rate vie the ILE-488 interface is 16 Ioes at 20 kHz kr a single channel, or 5 kHz for all three SQUID channels plus the auxiliary channel.

Analog Output: 4 analog outputs (60962) provided on the back panel for the 3 SQUID channels and the auxiliary analog input.

Autotuning of all SQUID parameters is cocoa, plished by single button push. All adushments may als be made manually or via the remote intoleces.


FLL Reset: Any channel may be reset manually or noternahcally et any user seksotable output voltage.

Bandwidth & Gain: Selectable bandwidths of 5 H. 509 HZ, 5 kHz & 50 kHz. (4-Pole Butterworth response). Selectable gains of (1, 2, 5, 10, ..., 500) corresponding to full-scale out.

### Flux Lock-Loop Waded IFL-301)

Two versions of the flux-locked loop are available. one for HIS sensors sod one for LTS sensors. The HIS FLL has a 25 kHz maximum bandwidth (selectable to be 250 Hz from the IMC-303 controller) and uses high-frequency bias reversal to minimize low-frequency noise intrinsic to the HTS sensors. This bias reversal eons not looresSe the white noise M the sensors or add any spikes within the specified bandwidth. The LTS FLL has a 53 kHz bandwidth (selectable to be 500 He from the IMC-333 controller) and uses no bias reversal since it is not required by the LTS sensors. Wider bendwidMs On both LTS and FITS leapt are available <sup>TM</sup> order.

All FLL functions are controlled remotely by the IMAG IMC-303 Controller. The FLU coenect to the IMC-393 via a cermpo- site cable. To minimize di, even when using vey lmo 011. bins, el high-frequency signals are transmitted by optical fiber between the FLL and Controller.



**TRISTAN TECHNOLOGIES** al AS Cosiereless  
Court East 5, Se 106, Can Grego C.4 82121

# SQUID iMAG

## DEWARs

- FEATURES
- Low scoldgef
- Close Tail Spacing
- 50 Standard and Custom Designs
- Comp
- lete
- Factory Testing
- Metall
- ic or
- Non-Metallic Construction

Tristen takes special pride in the Westale de. nen and constnsmn techniques n free developed. The use of SQUID magnetometers for tomegrolo or non-destructive lating and metie. (N010) measurements requires that magnetic signals front a subject al don temperature be coupled IC a superconducting pickup coil in the reets reserved II me dewar. It is easerreat to use nonmagnetic mandals and to have the smsted possible spacing between the oxygen reservoir and the outside of the dewar. Tularis development of adjustable tail covers have allowed marl was to be less than 2 me

PLMCI Series Liquid Helium InnumIrm  
Trietens IA113.10 is a fiberglass worm consigned Iv: magnetsm and Iff DE. The BMD-101,1 variant is supplied nor an upper aluminum bowing to resude weight, nonrotation CA, end Inredco reliability Intended for use with Tristan magnetometer probes. they provide a spacing of less than 10 mm be tween room temperature and ins kquid helium. The 6N610 typically uses 2 litreday of liquid helium. Custom deems with different size necks, tails, helium reservoirs or in-vacuum data:ellen coils are available.

NLO Series Liquid Nitrogen Dewar.


Speckaer anagnoned Ion use with HTS sbanu woers:sglean suers swine as-ortment of standard dewar designs. These include toted deware wely done access to the sensors , multi channel dewars, handheld deems pat eepe. In any ornantim and larger dewars ato more than 30-day hold does Chngemo inserts are noble to mount the SQUID caries rigidly in are clear and provide any performance features required of rho application eastern dewars wire mitered size necks, tails, or oxygen reservoirs can be spacial ordered.

### CONSTRUCTICIN TECHNIQUES AND MATERIALS

All "flown down incorporate loot egad daign cmrcnpts and are assembled cab the highest standards of workmanship. They are leak tested afar each phase of their construction end are cycled between 100 temperature and liquid name: Intorm:rakes macaws long-wen "liability. A complete series et tests is made at operating temperature including measurements of Ins equilibrium Outfall A lento, led report is supplied with each dewar.

The use of super-insulation end one or more vapornooled shields inlaid alimi news the need or liquid nitrogen in the BMD series. Tristarns own computer analysis is used to calculate the optimum layer deady of supreinspation In sect fern-arsine moon and the estaetan In careuhr applied by hand to main... Its density, even in thole Penal regions Ouch as corns, clow-spaced Cale or regfone where overtop occurs. Also computed are the number and position of do required ,ecor-coocied shields end, for custom steward the prodded cryogen boiled mi.

## DEWAR DIMENSIONS



### OPTIONS AND ACCESSORIES

- Dail-in-Vacuum dewars
- Adjustable ten option for liquid helium and NLD-500 series dewars. This option can alenral spacings kassthan 2 mm Iran liquid helium or re-heavn to room temperature
- tea Tilt Option for NLD Series 310 8 510 dewars
- Cumin designed low-pass (eddy current) filters end di shields.
- Mechanical anchoring of the helium reservoir to the outside Omar ease is available for applications where oachardal vibrations and relative mean may introduce noise
- Liquid helium or nitrogen level gauge. re Flexible metal transfer tubes wwwww awe did,1,7411,141,4



arm Cuweral ns Como East Sofa lee, San Easy\*. Ca 92121 Tv asS91 0601'00 Ass (SW 550-2 50 Ewer aSpereteseAS one



## do SQUID SENSORS

### Features:

- All Thin-Film Devices
- Niobium-Aluminum TH. Prifer Pi-ocean fm Robust LTS Device.
- YOCC Step-edge and thermostat Junctions fur Robust I-ITS Devices
- Symmetric Modulation Coils Eliminate Inductive Loading of Output

Tristan offers Revere! configurations of love...some SGI-110 sensors which serve as the heart of our iMAG SQUID systems-

Address your magnetic ...inning applications with the latest technology in both high-temperature and low.

The low-temperature 11TE7 SIMUILLiso run in liquid helium and are fabricated using a niobium/aluminum all thin-film hi-layer technology that combines durability with high sensitivity. Their Mature symmetric integral signal and modulation coils that eliminate output variations with varying input loads. The niobium-shielded package comes with screw terminals ready to accept your custom input circuit. Tristan can also provide thin-film integrated LTC SCUID magnetometers with slate-of-she.art performance. The Tristan Model L50125 can lee used with the Model SP Cryogenic Cable for ultrasensitive measurements of current (< (IT RANH.) and magnetic field to T fthHz). In conjunction with the Model [IMP and MFP Cryogenic Probes, it can measure a much wider range of electromagnetic properties in magnetic fields as high as 0 testa--see Tdstan's Cryogenic Probe den sheet for more Information.

The high-temperature (11-1T61 SOILPCLie run in liquid nitrogen at 77 K and are offered in magnetometer or gradiometer configurations. They feature Y13CD pick-up coils patterned on the chip and a tough passivation layer for protection from moisture and oxygen. All FITS MAO sensors use a common connector to attach them to the Model SP Cryogenic Cable; they may be easily interchanged to provide alternative pick-up coils and different sensitivity levels. We can guarantee magnetometer performance better than 90 ITOHz. For customers who need even knver noise levels and performance in magnetic fields we can pro-



### LTS Sensors

Input doll Inluoancla, TB pH

Noise level <5 a: 10<sup>-13</sup>/Hz



### HTS Sensors

Noise IHTM-11309: <90 fIAIHz


01 09-100ti e 10 pT/Hg

Sesides, the Standard LSO/90, Tristan can supply LTS sensors wven longer niobium shield cans such as that supplied with the MFP and RIMP probes (see Probe Sala Sheet for details). We can also supply the bare sensor chip for specialized applications.

The 2 pH input impedance 01 the LSCV20 allows eater matching of input circuit. Unlike asymmetric coil designs, the symmetric con design of the LSD/20 avoids inductive coupling of unwanted signals. Additionally, it exhibits no sensitivity variations with changing input impedances.

Tristing TITS sensors are the liron commercial devicaa to operate in both ambient and kilogauss environments. Step-edge junctions ensure uniform response Independent of sensor orientation, avoiding the Fraunhofer-like diffraction behavior seen in many monolithic bicrystal junction devices.

Tristan's FITS sensors are available in a wide variety of configurations. The standard FITS magnetometer sensor Is available in a 90° mounting (Model HTM-80) or Ina flexible end piece (Model HIM-1naba. The lee-tale section can be as long as 15 our without degrading performance. Pickup col dimensions other than the standard 5 mm X El ern are also available. The HTo-100 MiniMAG has a 50 in X 50 gilt pickup-coil and Is well suited for magnetic microscopy. The HTM-400's large 16 mm \* 15 min detection area gives it the highest sensitivity on any available HTS sensor. Tristan's gradiometers are available In either dElearr (shown be low) or cinvidg configurations.



Mama  
**M\*SO**

resolnann. (H\*)

**TRISTAN TECHNOLOGIES**

eves COMersterre Gond East Subt \*F. San cress, CA 92121  
reset asj:2tin pae: gum sso-zees 0,max: ma,enierero.ot

rees.enneen.hsse

## Ultra-high Resolution Scanning Magnetic Microscope

The Tristan model SMM-1000 Scanning Magnetic Microscope performs micron level non-destructive analysis of surface and sub-surface material properties using an array of small SQUID magnetometers. It can be used to Image diverse objects such as:

- micro-current distributions
- vortex motion in superconductors
- traces on a circuit board or multi-chip module
- weak electric currents in semiconductors
- integrated circuits
- magnetic domains



Model 1001 0100 Scanning SQUID Magnetometer

- The 65.1114-1000 us CH a proprietary intimated circuit That Incorporates an array of Superantducing Quantum Int-Heart.. Ideate. (SOLI IIts) to map the magnetic Cold from smell samples. 'rite use of liand helium SQUIDS provides a 100 cola improvement in sensitivity over other CIQIISPC detectors and allows high-resolution mapping of elect/to currents and magnetic sources located beneath the surface of the sample.
- It is a tally Natured measurement system that allows the user to extract a (napet le image of the object being measured Odor die entire dc - 10 kHz frequency range. Its flat phase response allows both in-phase and anadranne information to be oithined without distortion. Addliomd detection channels can Ire supplied to speed dela acquisition rates.
- It allows computer controlled earns ofobjects over alarge IS a 5 min) area with 0.17 Lcn stepping capability. The user has the ability to preprogram the scan coordinates.
- Automated setup and Computer control makes measurementa sapid and repeatable. System software provides the ability to control the critical system component, acquire data Horn the SQUID sensor, and analyze the data to date mine Me aingnaie proportion of the sample being measured. The use of open architecture software allows the user to modify and Customize nearly ol aspects or Operating including image accessisa.

**TRISTAN TECHNOLOGIES**

awpm. env daq on me ms Herechmone s. aa  
mums eef or .rno.laok mixing h 15  
rnoons-pourthey usco!

Besides measuring magnetic fields, the SMM-1000 can also be configured to detect:

- transient magnetic properties



SMM-1000 sample stage

### System Operation

The 01NM.1000 achieves micron resolution by the use of small ill use deamon coils and merge gag Fete.en rte coils ate the object(s) being camed.

The sample is mended inside an m change gas car et the lamer end of s cryogenic probe. This houses ell or the cryogenic portions of Me SW ant, lass) a measurement, Is Sled with a man amount of he/enges. The sample to placed on Mesmeric stage and the erotic can ettachoel. Then ITreeMO.1 Prete Is lowered mu He fined helium dewy Stoe the sample arose lmo armed IO \*2



otos er sq. tower fliv onusescene re ren Selectee coil\*

### System Component's

The standard model SM51.1000 pldoes a single channel WM probe nlagnete Ontonen autsertem) end WAG\* 891013 e..11311.. UMW p5111 measurement end control gobyyet,

Well hem, gear veil vitromen stand, onoto wenninis, 31.1 gas dandling states monomer cont. console and complete software package for system control, des acquisition end dele anlimes me model 7011.4111110 on he supplied and additional. capettillos ...extend Its measurement capabilities.

### OPTIONS AND ACCESSORIES

**Addionie Defective chgeortekt** The model ,SPAM-101774 mensuremeal opethites can be untended to wfil-channel twprahnev Adratonel mire. (Bit nenswomell sites 050 ts3 installed le reduce massuremo. ems The stereo, stance beams, the LO118 is 54 gm. coils maybe larded 10C pm, 120 rim, or 209 era span at no este charge Other .41 diameters and nonNitalurions are ayesshoe as options.

**Spilled Feld** This option generates a vertical dri, a dc

magnetic MI on the sample This snows mgngrel encsplicity measuremerle on insulators, conductors anal ferrous materials to he deleared.

**Variable swept\* rerper.111:** The Wandard measurement temperatrise is 42 IC The variatla lemoacature \*peon elates sample temperatrise to Lis vated between 2K and los


### SPECIFICATIONS

**SENSOR:** Lore tbleEBMLNP weemomshring SRIFJfke, interference device (512110)Ortratea .92X

**SPATIAL RESOLUTION:** I Onto single spate sources

**SENSITAITro** leSatits [WC pT.Hr)

**OPE/1617HG NANO/ADM de** - 10 kHz. M0.1.11.10nrs mn bar made at any frequency.Bandeoidths above to kHz ere aveffleble



**TRISTAN TECHNOLOGIES**