



A new, noninvasive investigational tool for pre- and full-term infants



- Mapping of sites and dynamics of sensory functions - auditory, somatosensory, and visual modalities
- Assay stages of nervous system development



Somatic evoked magnetic field (SEF) obtained from a 7-month old as a function of number of averages from N=4 to 173 epochs. The waveforms are the differences of the SEF at two field extrema. This shows that a small number of averages are needed to acquire SEF data. (data acquired using a Tristan babySQUID[®] system).

Unique Features of Artemis123[®]

Superior spatial resolution and sensitivity

- Significantly more sensitive to neuronal sources than conventional whole-head MEG systems
- Similar or better spatial resolution compared to existing whole-head MEG sensors
- Better spatial resolution than EEG (EEG signals are distorted by skull defects (fontanels and sutures), making it difficult to localize epileptiform tissue)
- Rapid scanning: a typical clinical scan can be completed within thirty minutes
- Anti-vibration construction; infant motion will not cause vibrational artifacts
- Sensor noise level $< 10 \text{ fT}/\sqrt{\text{Hz}}$
- A dense array of closely-spaced sensors located just below the outer surface of a headrest.
- Allows simultaneous measurement of the occipital area and parietal and temporal areas
- Includes position tracking device and software, permitting measurements during sleep or relatively quiescent wakefulness
- The measurement cradle and companion electronics cart are portable and can be wheeled in and out of elevators, obstetric suites and neonate ICUs

Mapping brain function and detecting neurological abnormalities in infants

Detection of cortical function in newborns is needed for clinical intervention in the early stages of neurological disorders, before external signs appear and the conditions develop and worsen. Areas where Artemis123^{®1} can be used for neonatal neurological assessment include:

- Autism
- Epilepsy
- Cerebral palsy
- Perinatal asphyxia
- Hypoxic-ischemic encephalopathy
- Periventricular white matter injury
- Monitoring recovery from trauma

Identifying how infants learn is of interest to many sectors of society, but such studies rely heavily on behavioral analyses. Having a direct measure of cortical activity would provide precise information on the dynamic response in the brain during learning processes. Potential uses of Artemis123[®] for developmental studies include:


¹ Artemis123[®] is a registered trademark of Tristan Technologies, Inc. All rights reserved



Artemis123[®] System Description

Principles of Operation

Like adult Magnetoencephalography (MEG) systems, Artemis123[®] uses superconducting sensors to non-invasively detect and map magnetic fields generated by cortical neural activity. However, Artemis123[®] takes advantage of the fact that the infant's scalp and skull are very thin. Tristan's fabrication methods put the sensing coils very close to the infant brain's sources of activity, even though SQUIDs must operate in an ultra-cold liquid helium environment. The net result is a significant increase in amplitude of neonate MEG signals. Also, the high density of detectors results in higher spatial resolution compared to adult whole-head MEG.




System Components

- Sensor/Cradle/Bed on mobile cart – easily accessed height
- Power supplies and computer on companion mobile cart to minimize noise
- Subject Tracking - optical tracking system updates movement at 30 Hz with ½ mm accuracy
- Part-wise mapping or optional optical one-click 3D imaging system

SQUID Sensor Array

- 606 cm² sensor coverage area
- 100+ detection coils
- Coil type: 15 mm-diameter first order gradiometers. Adjacent coils can be electronically combined to form planar gradiometers
- Coil gap: ~8 mm from sensor to outer surface
- Coil sensitivity: better than 10 fT/√Hz
- Reference channels: 12-element tensor array for noise reduction

BabySQUID[®] Patient Data



A 9-year old with cerebral palsy (unilateral congenital malformation in the left hemisphere in the lower somatosensory and motor cortices). Left: MRI image of the malformation. Right: Somatosensory evoked magnetic field (SEF) produced by tactile stimulation to the tip of left and right index fingers. Note a strong SEF in the intact right cortex after left index finger stimulation, and abnormal small initial cortical response and ipsilateral response in the affected left hemisphere after the right index stimulation.

Power and Physical Requirements

- Power: 1.5 kW filtered circuit
- Patient bed: 1 m x 2 m x 1.1m (40" x 79" x 42")
- Patient bed weight: 200 kg (440 lbs)
- Instrument cart weight: 150 kg (330 lbs)

Larger coverage areas, higher channel counts, and/or different coil dimensions and configurations are available on request. Contact Tristan for additional information.

All Tristan products are covered by a 1-year warranty. Service contracts may be purchased to provide post-warranty coverage.



6185 Cornerstone Court East, Suite 106
 San Diego, CA 92121
 (858) 550-2700 [fax] 550-2799
<http://www.tristantech.com>

The Tristan babySQUID[®] (patents issued and pending) and Artemis123[®] are classified as investigational devices and are currently offered for research use only. Tristan is in the process of seeking both CE (European) medical device directive and FDA (U.S.) certification for clinical use. Specifications are subject to change without notice.